

## MS-PA

### LED Digital Magnetic Stirrer

With excellent stirring functionality, the compact shape of FlatSpin takes less bench space and makes storage easier.

## MS-PB

### Magnetic Stirrer

## Features

- Digital speed control within the range of 0–1500rpm
- Maximum stirring quantity of H<sub>2</sub>O at 3L
- LED display shows stirring speed (MS-PA)



Specifications	MS-PA	MS-PB
Work plate Dimension	Ø135mm ( 5 inch )	Ø135mm ( 5 inch )
Work plate material	Nylon+GF	Nylon+GF
Motor type	Brushless DC motor	Brushless DC motor
Stirring positions	1	1
Max. stirring quantity[H <sub>2</sub> O]	3L	3L
Max.magnetic bar[length]	50mm	50mm
Speed range	0–1500rpm	0–1500rpm
Speed display	LED	scale
Protection class	IP42	IP42
Motor input power	5W	5W
Motor output Power	3W	3W
Power	15W	10W
Voltage	100–120/200–240V, 50/60Hz	100–120/200–240V, 50/60Hz
Dimension[WxDxH]	150x260x80mm	150x260x80mm
Weight	1.8kg	1.8kg
Permissible ambient temperature and humidity	10–40°C,80%RH	5–40°C,80%RH



Strong magnetism



Chemical resistance



LED display